Činidlo  
pro test ACETOIN

Kat. č.: MLT00016

Pro mikrobiologii

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Cat. No.: MLT00016

For microbiology

The Reagent for ACETOIN test is supplementary colour forming reagent for the test for the formation of acetoín (the Voges-Proskauer reaction) included in some of MIKROLATEST® kits. The instructions for the use of reagent are a part of the working instructions of the respective MIKROLATEST® kit.

**Preparation for use:**

Pour the content of the bottle labelled Reagent for VPT I – solvent into the bottle labelled Reagent for VPT I – weight (α-naphtol). After the content has dissolved, pour carefully the solution arisen into the dropping bottle labelled Reagent for VPT I. It is necessary to record immediately the DATE OF PREPARATION and EXPIRY (see Storage, expiration) in the bottle label.

The dropping bottle labelled Reagent for VPT II is ready to use.

**Reagent for ACETOIN test contains:**

Reagent for VPT I – weight	1,8 g
Reagent for VPT I – solvent	18 ml
Reagent for VPT II	18 ml
Dropping bottle for Reagent for VPT I	1 pc
Information leaflet	

**Storage, expiration:**

Store the package of Reagent for the ACETOIN test in the dark at the temperature of (+2 to +25) °C. Expiration is marked at each package.

The dropping bottle with the solution for VPT I must be stored in the dark at the temperature of (+2 to +8) °C, the prepared solution is usable up to 4 months since the date of preparation – not to exceed the expiry date of the kit.

The dropping bottle with reagent for VPT II can be stored at the temperature of (+2 to +25) °C.

**Warning:**

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Reagents of the kit contain ethanol (95 %, VPT I – solvent), methanol (5 %, VPT I – solvent), potassium hydroxide (45 %, VPT II), and α-naphtol (100 %, VPT I – weight).

**Reagent for VPT I – weight**

H302+H312 Harmful if swallowed or in contact with skin.

H315 Causes skin irritation.

H318 Causes serious eye damage.

H335 May cause respiratory irritation.

P261 Avoid breathing vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P313 IF exposed or concerned: Get medical advice.

Danger

**Reagent for VPT I – solvent**

H225 Highly flammable liquid and vapour.

H302 Harmful if swallowed.

H371 May cause damage to organs (eyes).

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P261 Avoid breathing vapours.

P280 Wear protective gloves/protective clothing/eye protection.

P308+P313 IF exposed or concerned: Get medical advice.

Danger

**Reagent for VPT I (dropping vial)**

H225 Highly flammable liquid and vapour.

H302 Harmful if swallowed.

H318 Causes serious eye damage.

H371 May cause damage to organs (eyes).

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P261 Avoid breathing vapours.

P280 Wear protective gloves/protective clothing/eye protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**Reagent for VPT II**

H302 Harmful if swallowed.

H314 Causes severe skin burns and eye damage.

P280 Wear protective gloves/protective clothing/eye protection.

P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P313 IF exposed or concerned: Get medical advice.

Danger

**Waste disposal:**

Contingent rest of the reagent should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and containers to recycling.

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