Kat. č.: MLT00002

**Pro mikrobiologii**

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Cat. No.: MLT00002

**For microbiology**

The kit CANDIDAtest 21 is designed for the routine identification of the most clinically relevant yeasts. The kit enables identification of thirty-four strains by means of twenty-one biochemical tests (chromogenic substrates, decarboxylase and assimilations). The tests are placed in wells of microtitration plates. Three rows, each consisting of eight wells, are intended for the identification of one strain. Negative tests for each group of reactions are placed on the plate to help with the correct readings of positive results.

**CANDIDAtest 21 kit contains:**

- 5 microwell plates with tests (each for the identification of 4 isolates)
- Instructions for use
- Colour scale for CANDIDAtest 21
- 20 record sheets
- Lid
- 20 foils for incubation
- Storage bag for an open plate
- 20 suspension media in glass test tubes

**Storage, expiration:**

- The CANDIDAtest 21 kit and the suspension media should be stored at (+15 to + 25)° C in original packaging. Date of expiry is indicated on each package. After removal of a microplate from its packing, use the microplate within 4 weeks.

**Recommended instructions for use for CANDIDAtest 21**

**Material required to perform a test (not included in the kit):**

- Suspension medium for CANDIDAtest 21
- Petri dishes with the cultivation medium (Sabouraud-dextrose agar 2% without additives)
- Instrument DENSILAMETER II (Cat. No. INS00062)
- Automatic micropipette 0.1 ml, sterile tips
- Thermostat 30°C
- Common microbiological laboratory equipment (loops, markers, burner)

**For results evaluation (not included in the kit):**

- Code Book for CANDIDAtest 21 - located at [www.erbalachema.com](http://www.erbalachema.com)
- The ErbaExpert Identification Program

**Caution:**

For professional use only

**Respect the rules for work with infectious material!**

**Isolation of cultures:**

- The isolation of culture should be carried out by usual techniques on Sabouraud-Dextrose-(2%) agar (without additives)

**Preparation of inoculum:**

- Before you start working, rock the content of suspension medium for CANDIDAtest 21 gently. Using an 24 h culture (48 h for slowly growing species) prepare a microbial suspension in the suspension medium. The suspension must have turbidity equal to McFarland 0.5. Lower or higher turbidity can cause false results. Avoid quick movements during the preparation of suspension which can cause formation of air bubbles and thus influence the measurements of suspension density. Use rocking movements or a loop to homogenise microbial suspension.

**Purity control:**

- Confirm purity of the suspension by inoculation of the yeast suspension on a Sabouraud-Dextrose-(2%) agar plate and incubate for 24 hours. Check the purity of the microbial growth. Incubate plates with weak growth for further 24 hours.

**Preparation of CANDIDAtest 21 plate:**

- Cut off one side of the packaging foil and remove the plate from the package.
- Cut off the required number of the strips from the plate (3x8 wells per identification) and insert them into a prepared empty frame. If you work with a MIKRO-LA-TEST for the first time and no empty frame is available, use the frame of the first plate. Insert any unused strips into a storage bag freely.
- Write identification numbers of tested strains on the corresponding strips.

**Note:**

Do not place individual strips too close to each other to avoid contamination of adjacent wells by inoculated suspension.

- Insert remaining strips of microplate to the storage bag and keep it in a dry place at room temperature. Don't store it longer than for 4 weeks.
- Clean the inner side of the lid by ethanol before you use it to cover a plate during your work. Disinfect the frame and the lid if you wish to re-use it.

**Note:**

Any uneven distribution of substrate in the well does not affect the functionality of the test.

**Inoculation:**

- Inoculate 0.1 ml of the suspension into all wells of the respective three row-strip.

**Incubation:**

- Cover the test strips with the incubation foil.
- Incubate CANDIDAtest 21 plate at (+25 to +30)°C for 24 hours and for 48 hours for species growing slowly.

**Reading:**

- Evaluate reactions after 24 hours or 48 hours of the incubation.
- Remove the incubation foil if necessary.
- Read all the results and record the results in the report sheets. Read the assimilation reactions against black text or a black line on a white background.

Note: Read the reactions in accordance with the Colour scale for CANDIDAtest 21, the table "Interpretation of reactions" and/or according to the colour reactions of control strains.

**Identification:**

- Calculate a 7-digit numerical code: The tests are divided into the groups of three. Negative reactions are always assigned the numerical value 0. Positive reactions are assigned a value as marked in the report sheet. The identification is performed by summing up all the positive test values in the group and creating a 7-digit number. Find the code number in the Electronic Code Book.

- All the taxa used in the Electronic Code Book point to the imperfect form of a respective yeast. Please note that the use of taxa nomenclature of perfect form is only justified when sexual spores are proven. Physiological tests cannot distinguish between sexual and asexual form.

**Additional tests:**

Remarks may appear below the test results in the Electronic Code Book - in order to indicate an insufficient growth or to specify additional discriminatory tests.

**GET – Germ tubes test**

- Inoculate yeast colonies in 1.0 ml of human or animal serum (density 0.5-1 McF).
- Place it in an incubator at (+35 to +37)°C for 3 hours.
- Observe germ tube formation under a microscope (magnification: 100 x): Germ tubes of genus *Candida* grow in a filamentous form. The formation of a filamentous form of germ tubes is not specific for the genus *Candida* (eg., *C. albicans*, *C. dubliniensis*, *C. africana*) after an overnight incubation.

**PSH – Pseudohyphes**

A characteristic of *Candida* is the detection of pseudomycelium, which will be formed in the presence of nutritionally poor substrates.

- Prepare rice extract agar according to instructions for use of a supplier or a standard microbiology protocol.
- Inoculate a yeast strain onto a thin rice agar plate and cover with a coverslip.
- Place it in an incubator at (+25 to +30)°C for 4 days. Observe biomass growth under a microscope (magnification: 100x) every day.

**Evaluation:** Pseudomycelium consisting of pseudohyphes with terminal chlamydo spores (spores with a spherical shape and a refractive cell wall) is indicative for *Candida albicans* and *C. dubliniensis*. Most of the other *Candida* sp. form pseudomycelium consisting of pseudohyphes without chlamydo spores. However some species do not form pseudohyphes.

**ATS – Arthrospores**

A characteristic of genus *Geotrichum* and *Trichosporon* is formation of arthrospores. Arthrospores are formed by the fragmentation of the terminal parts of hyphae.

**Disposal of used material:**

- All ampoules, tips and strips must be autoclaved or incinerated after the use.
- Discard paper packaging waste to recycling.

**The most frequent causes of identification failure:**

- Contaminated culture: work with pure cultures isolated from Sabouraud-Dextrose-(2%) agar (without additives) not older than 24 hours, 48 hours for slowly growing yeast. The Sabouraud-Dextrose agar was used for database setup. Usage of different culture media can affect the reaction profile and therefore interfere with the phenotypic identification.
- Unsuitable storage conditions: Store at temperature (+15 to +25)°C.
- Used inoculum was of low density or small volume: please follow the correct Mc Farland 0.5 density of the suspension. Pay attention to the homogeneity of the inoculum.
- Inoculum has contaminated adjacent strips.
- Cover the strips only with the manufacturer's sealer.
- Follow the exact incubation time of 24-48 hours.
- Failure to observe some of the steps in the recommended procedure.
- There may be a species or a strain not included in the Taxa list.

Note: A pouch filled with the indicator "silica gel" is supplied with the sealed test plates. A colour change of the indicator pouch from blue to pink may indicate trapped moisture. Please do not use this plate.

**Performance:**

The kit was tested on a set of 218 strains and the results were compared with a reference method. The identification results of both methods corresponded for 192 strains. 9 strains were not identified.

**Quality control of CANDIDAtest 21:**

The plates and reagents are systematically controlled for their quality at various stages of the production. The manufactured batches are tested using standard bacterial cultures. For those who wish to perform their own quality control tests, the following cultures are recommended:

- *Candida albicans* CCM 8261 (ATCC 90028)
- *Cryptococcus neoformans* CCM 8312 (ATCC 90112)

The strains are supplied in freeze-dried ampoules by the CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz).

**Caution:**

It is necessary to use fresh isolates of the CCM strains each time when a kit is tested for its functionality. **These strains are recommended to check the functionality of the kit and not to monitor accuracy or effect of the identification.**

Row	H	G	F	E	D	C	B	A
<b>Candida albicans CCM 8261</b>								
1	NGA	αGLU	βGLU	ONC	GENA	GLU	GAL	MALA
	+	+	+		–	+	+	+
2	αGAL	PHE	URE	UCO	RHA	INO	TRE	LAC
	–	v	–		–	–	+	–
3	PRO	PGUR	MEL	ACO	XYL	CEL	SUC	RAF
	+	–	–		+	–	+	–
<b>Cryptococcus neoformans CCM 8312</b>								
1	NGA	αGLU	βGLU	ONC	GENA	GLU	GAL	MALA
	–	+	+		+	+	+	+
2	αGAL	PHE	URE	UCO	RHA	INO	TRE	LAC
	–	+	+		+	+	+	–
3	PRO	PGUR	MEL	ACO	XYL	CEL	SUC	RAF
	–	+	–		+	s	+	s

Explanations: + = Positive reaction    – = Negative reaction    s = Weak reaction  
v = Variable reaction

**CANDIDAtest 21 INTERPRETATION OF REACTIONS**

Column	Test	Code	Reaction	
			positive	negative
Row 1				
H	N-Acetyl-β-D-Galactosaminidase	NGA	yellow	colourless
G	α-Glukosidase	αGLU	yellow	colourless
F	β-Glukosidase	βGLU	yellow	colourless
E	Control of chromogenic substrates	ONC	Negative control	
D	Gentiobiose	GENA	growth	no growth
C	D-Glucose	GLU	growth	no growth
B	Galactose	GAL	growth	no growth
A	Maltose	MALA	growth	no growth
Row 2				
H	α-Galactosidase	αGAL	yellow	colourless
G	L-Phenylalanine-aminopeptidase	PHE	yellow	colourless
F	Urease	URE	red, red-orange	yellow, yellow-orange
E	Urease control	UCO	Negative control	
D	L-Rhamnose	RHA	growth	no growth
C	Inositol	INO	growth	no growth
B	Trehalose	TRE	growth	no growth
A	Lactose	LAC	growth	no growth
Row 3				
H	L-Prolinaminopeptidase	PRO	yellow	colourless
G	p-Nitrophenyl-β-Glukuronidase	PGUR	yellow	colourless
F	Melibiose	MEL	growth	no growth
E	Assimilations control	ACO	Negative control	
D	D-Xylose	XYL	growth	no growth
C	Cellobiose	CEL	growth	no growth
B	Saccharose	SUC	growth	no growth
A	Raffinose	RAF	growth	no growth

**Limitation:**

The kit CANDIDAtest 21 is designed only for identification of the species included in the database (see Taxa list). Other microorganisms can not be identified by the kit.

**Taxa List**

▪ <i>Candida africana</i>	▪ <i>Candida krusei</i>	▪ <i>Candida valida</i>
▪ <i>Candida albicans</i>	▪ <i>Candida lambica</i>	▪ <i>Cryptococcus albidus</i>
▪ <i>Candida catenulata</i>	▪ <i>Candida lipolytica</i>	▪ <i>Cryptococcus humicola</i> Komplex
▪ <i>Candida dubliniensis</i>	▪ <i>Candida lusitaniae</i>	▪ <i>Cryptococcus neoformans</i>
▪ <i>Candida famata I</i>	▪ <i>Candida pelliculosa</i>	▪ <i>Cryptococcus terreus</i>
▪ <i>Candida famata II</i>	▪ <i>Candida membranefaciens</i>	▪ <i>Geotrichum candidum</i>
▪ <i>Candida famata III</i>	▪ <i>Candida norvegensis</i>	▪ <i>Geotrichum capitatum</i>
▪ <i>Candida famata IV</i>	▪ <i>Candida norvegica</i>	▪ <i>Rhodotorula glutinis</i>
▪ <i>Candida glabrata</i>	▪ <i>Candida parapsilosis</i>	▪ <i>Rhodotorula mucilaginosa</i>
▪ <i>Candida guilliermondii</i>	▪ <i>Candida magnoliae</i>	▪ <i>Saccharomyces cerevisiae</i> TRE -
▪ <i>Candida inconspicua</i>	▪ <i>Candida rugosa/pararugosa</i>	▪ <i>Saccharomyces cerevisiae</i> TRE +
▪ <i>Candida intermedia</i>	▪ <i>Candida tropicalis</i>	▪ <i>Trichosporon species</i>
▪ <i>Candida kefyr</i>	▪ <i>Candida utilis</i>	▪ <i>Trichosporon species</i> RAF - /MEL -

**Health protection:**

Components of the kit are not classified as dangerous.

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