Kat. č.: MLT00003

Pro mikrobiologii

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Cat. No.: MLT00003

For microbiology

The diagnostic kit CANDIDA-Screen is designed to screen the most frequent clinically significant yeast species. CANDIDA-Screen is placed in wells of a one-strip microtitration plate. It is possible to carry out 12 identification on a plate.

**CANDIDA-Screen kit contains:**

- One-strip microtitration plates with desiccant, 3 pc
- Instructions
- Three PE bags for incubation
- A storage bag for unused plate
- Colour scale for CANDIDA-Screen kit
- Record form
- Lid

**Storage, expiration:**

CANDIDA-Screen should be stored at (+2 to+8)° C. Date of expiry is indicated on each package.

**Safety rules:**

The CANDIDA-Screen kit is for professional use only. Only a trained person should work with the kit. This person should follow safety rules for work with infectious material and for safe disposal according to the guidelines of the workplace.

### Recommended instructions for use of CANDIDA-Screen

**Material required to perform a test (not included in the kit):**

- Petri dishes with the cultivation medium Sabouraud-dextrose agar 2% without additives, Columbia blood agar or another suitable non-selective cultivation medium.
- Paraffin oil sterilised (Cat. No. MLT00042 - 40 determinations)
- Tubes (100x15 mm) with sterile unbuffered physiological solution
- Instrument DENSILAMETER II, cat. No. INS00062
- Vortex V1( Cat. No. 50001715)
- Multistep pipette Mikrolastepper, (Cat. No. 50001707)
- Incubator 25 - 30 °C
- Common microbiological laboratory equipment (loops, burners, markers)

**For results evaluation (not included in the kit):**

- Code Book for CANDIDA-Screen - located at [www.erbalachema.com](http://www.erbalachema.com)
- The ErbaExpert Identification Program

**Isolation of cultures:**

- The isolation of culture should be carried out by usual techniques on Sabouraud-Dextrose 2% agar (without additives) or blood agar. Another suitable culture medium can be used after it was tested on recommended control strains.

**Preparation of inoculum:**

Prepare microbial suspension from pure well grown 24 hour (48 hours for slowly growing yeast) culture in unbuffered sterile physiological solution. The suspension must have turbidity equal to McFarland 3. Use vortex to homogenise the suspension well.

**Purity control:**

- Confirm purity of the suspension by inoculation of the yeast suspension on Sabouraud-Dextrose agar (2%) using the same loop as in the preparation of the inoculum. Incubate it for 24 hours. Check the purity of the culture after 24 hours. Incubate for another 24 hours if the growth is weak.

**Preparation of CANDIDA-Screen plate:**

- Cut off one side of aluminium packaging foil and remove the plate from the package. Separate required number of strips of microtitration plate. Remove the adhesive tape from individual strips and insert them into prepared frame. In case you work with MIKROLATEST® kit for the first time and an empty frame is not available, use the frame of the first plate. The unused strips of the first plate put into the storage bag freely.
- Insert any unused strips into enclosed aluminium bag with desiccant and store it at the temperature (+2 to +8)° C for further use. Protect it against the moisture. We recommend to use the plate before the 4th week of its opening.

**Note:** Any uneven distribution of substrate in the well does not affect the functionality of the test.

**Inoculation:**

- Inoculate 100 µl of the suspension in the physiological solution into the wells of a strip.
- Add 3-4 drops of paraffin oil into inoculated wells.

**Incubation:**

- Insert the inoculated microtitration plate into an incubation PE bag and fold the open end of the bag under the plate to prevent drying of the inoculum.
- Incubate the Candida-Screen plate in an incubator set up for 25 - 30°C for 24 hours.
- It is possible to extend the incubation for slowly growing organisms such as *C. neoformans* for further 24 hours.

**Reading:**

- Evaluate reactions after 24 hours (eventually after 48 hours).
- Remove the plate from the incubation bag if necessary.
- Read all the reactions and record the results in the report sheets.

**Note:** Read the reactions according to the Colour scale or the table **Interpretation of the reactions** or according to the colour reactions of the control strains. For evaluation use **Code book** or use software ErbaExpert.

**Additional tests:**

The identification of the yeasts according to the **Differentiation table** can give indefinite results in some cases. We recommend to perform the following additional tests:

**GET - Germ tubes test**

- Inoculate yeast colonies in 1.0 ml of human or animal serum (density 0.5-1 McF)
- Place it in an incubator at 35–37 °C for 3 hours.
- Observe formation of germ tubes under microscope (magnification 100x). Germ tubes of genus *Candida* grow in filamentous form. The formation of filamentous form of germ tubes is not specific for the genus *Candida* after an overnight incubation. A positive result of the test for the species included in the **Differentiation table** applies to *C. albicans* only

**PSH - Pseudohyphae**

A characteristic for genus *Candida* is detection of pseudomycelium which will be formed in the presence of nutritionally poor substrates.

- Prepare rice extract agar according to the instructions for use of a supplier or a standard microbiology protocol.
- Inoculate a yeast strain onto a thin rice agar plate and cover with a cover slip.
- Place it in an incubator at 25–30 °C for 4 days. Observe a biomass growth under a microscope (magnification 100x) every day.

**Evaluation:** Pseudomycelium consisting of pseudohyphae with terminal chlamydo spores (spores with a spherical shape and a refractive cell wall) is indicative for *Candida albicans*. Most of other *Candida* sp. form pseudomycelium consisting of pseudohyphae without chlamydo spores. However some species don't form pseudohyphae.

**HYP - Formation of true hyphae**

A characteristic of genus *Geotrichum* and *Trichosporon* is formation of true hyphae and arthrospores. Arthrospores are formed by the fragmentation of the terminal parts of the hyphae.

## Interpretation of the reactions:

Column	Test	Abbreviation	Reaction: positive	negative
H	Urease	URE	red, red-to-orange	yellow, light orange
G	Sucrose	SUC	yellow, yellow-to-green	green
F	Maltose	MLT	yellow, yellow-to-green	green
E	Lactose	LAC	yellow, yellow-to-green	green
D	Galactose	GAL	yellow, yellow-to-green	green
C	Trehalose	TRE	yellow, yellow-to-green	green
B	Celobiose	CEL	yellow, yellow-to-green	green
A	Proline	PRO	yellow, yellowish	colourless

### Kit characteristics:

The kit was tested on a set of 60 strains. The CANDIDA-Screen identified 50 strains correctly (83.3%).

### The most frequent causes of identification failure:

- Contaminated culture: work with pure cultures isolated from recommended culture media. Culture should be 24 hours old (48 hours for slowly growing species).
- Used inoculum was of low density or small volume. Please follow the correct McFarland 3 density of the suspension. Pay attention to the homogeneity of the inoculum.
- Inoculum has contaminated adjacent strips.
- Failure to observe some of the steps in the recommended procedure.
- There may be an atypical strain or a species not included in the Taxa list.

### Disposal of used material:

Insert used strips into a container for infectious material and autoclave it or incinerate it.

### Quality control of CANDIDA-Screen kit:

The following strains are recommended for internal quality control:

**Candida albicans CCM 8261 (ATCC 90028)**

**Cryptococcus neoformans CCM 8312 (ATCC 90112)**

**Klebsiella pneumoniae subsp. pneumoniae CCM 5852 (ATCC 13882)**

These strains are supplied by the CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz). The strains are supplied lyophilised or on gelatin discs.

### Caution:

It is necessary to use fresh isolates of the CCM strains each time when a kit is tested for its functionality. These strains are recommended to check the functionality of the kit and not to monitor accuracy or effect of the identification.

### Control strains:

Row	H	G	F	E	D	C	B	A
<b>Candida albicans CCM 8261</b>								
1	URE	SUC	MLT	LAC	GAL	TRE	CEL	PRO
	-	+	+	-	+	v	-	+
<b>Cryptococcus neoformans CCM 8312</b>								
1	URE	SUC	MLT	LAC	GAL	TRE	CEL	PRO
	+	-	-	-	-	-	-	-
<b>Klebsiella pneumoniae CCM 5852</b>								
1	URE	SUC	MLT	LAC	GAL	TRE	CEL	PRO
	-	+	+	+	+	+	+	-

### Explanations:

+ = positive reaction

- = negative reaction

v = variable reaction

### Differentiation table:

Taxon	URE	SUC	MLT	LAC	GAL	TRE	CEL	PRO
<b>Candida albicans</b>	-	+	+	-	+	d	-	+
<b>Candida glabrata</b>	-	-	-	-	-	+	-	-
<b>Candida guilliermondii</b>	-	+	-	-	+	v	(-)	+
<b>Candida kefyr</b>	-	+	-	(+)	+	-	v	-
<b>Candida krusei</b>	-	-	-	-	(-)	-	-	(-)
<b>Candida lipolytica</b>	(+)	-	-	-	(-)	-	-	+
<b>Candida lusitanae</b>	-	+	(+)	-	v	v	+	+
<b>Candida parapsilosis</b>	-	+	v	-	+	v	-	+
<b>Candida tropicalis</b>	-	+	+	-	+	+	-	-
<b>Cryptococcus neoformans</b>	+	v	(-)	-	v	(-)	(-)	-
<b>Geotrichum sp.</b>	-	-	-	-	v	-	-	(-)
<b>Saccharomyces cerevisiae</b>	-	+	(+)	-	v	v	-	-
<b>Trichosporon sp.</b>	(+)	-	-	-	-	-	-	(-)

Explanations: + = positive reaction - = negative reaction v = variable reaction (+) = mostly positive reaction (-) = mostly negative reaction

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