Kat. č.: MLT00010

Pro mikrobiologii

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Cat. No.: MLT00010

For microbiology

The NEFERMtest 24 kit is designed for routine identification of Gramnegative non-fermenting bacteria and families *Vibrionaceae*, *Aeromonadaceae* and *Pleiomonas shigelloides* present mainly in clinical material. The kit enables identification of forty strains using twenty four biochemical tests. The tests are placed in the wells of a segmented microtitration plate. A strain is identified by tests placed in three rows with eight wells each. The Oxitest is an additional („off-line“) test that is essential for the identification.

A test for fermentation of glucose, supplied as OFtest in microtitration wells, can be used to confirm that the isolate belongs to non-fermenting bacteria.

The NEFERMtest 24 kit contains:

- 10 microtitration plates (each for identification of 4 strains) with desiccant
- Instructions for use with a differentiation table
- 40 record sheets
- 10 PE bags for incubation
- A storage bag (to insert any unused strips), 1 pc
- Lid
- Colour scale for NEFERMtest 24 kit

Storage, expiration:

NEFERMtest 24 must be stored at the temperature (+2 to +8) °C. The date of expiration is indicated on each package.

Recommended instructions for use for NEFERMtest 24

Material required to perform NEFERMtest 24 that is not included in the kit:

- Paraffin oil sterilised (Cat. No. MLT00042, 150 determinations)
- Petri dishes with cultivation medium (blood agar)
- Test tubes with 3 ml of sterile **unbuffered** physiological solution
- Instrument DENSILAMETER II (Cat. No. INS00062)
- Automatic micropipette 0.1 ml, sterile tips
- Incubator 37 °C
- Regular microbiological laboratory equipment (loops, markers, burner)

Recommended test to confirm that an isolate belongs to non-fermenting bacteria:

- OFtest (Cat. No. MLT00032, 288 determinations)

Additional („off-line“) tests :

- Oxitest (Cat. No. MLT00039, 50 determinations)
- Reagent for OXIDASE test (Cat. No. MLT00022, 250 determinations: it highlights the reaction)

Note: It is necessary to use Oxitest for correct identification. This test is necessary to create a profile for identification by the Code book as well.

Identification aids:

- Code Book for NEFERMtest 24 - located at www.erbalachema.com
- The ErbaExpert Identification Program

Caution:

- For professional use only

Respect the rules for work with infectious material!

We recommend to verify colour expression of positive and negative reactions on standard control organisms (see below) if you work with the kit for the first time. Further, we recommend to verify the colour of the reactions with every new batch of NEFERMtest 24.

Isolation of cultures:

- Isolation of the cultures is carried out by conventional bacteriology methods on the blood agar. To confirm that the isolates belongs to non-fermenting bacteria, perform a test for glucose fermentation OFtest (the test is read in max. 4 hours)

Preparation of inoculum:

- Prepare a bacterial suspension in a 3 ml unbuffered saline solution from a pure 24 hour culture (or 48 hour culture for slowly growing bacteria). Homogenise the suspension well.
- The suspension must have turbidity equal to the grade 2.0-2,5 of McFarland turbidity scale. The suspension of lower or higher turbidity can cause false reactions.

Performing additional test Oxitest:

- Perform an Oxitest from a 24 hour culture and record the results into a record sheet.

Confirmation of inoculum purity:

- To confirm the purity of the tested inoculum, streak out a sample from the prepared suspension using the same loop as in its preparation process. Check the culture purity after 24 hours of incubation. Incubate the plates with a weak growth for another 24 hours. Purity control culture can be used to perform additional tests as well.

Preparation of NEFERMtest 24 plate:

- Open an aluminium bag by cutting it close to the weld and take out a plate.
- Cut off a required number of rows (strips) of a plate using a scalpel in accordance with a number of tested strains (three rows, i.e. 3x8 tests for identification of one strain).
- Remove the strips from the panel, remove the adhesive aluminium foil from the strip and insert it into an empty frame. If you are working with MIKRO-LA-TEST® for the first time and an empty frame is not available, use the frame of the first plate. Insert any unused strips of the first plate into a storage bag freely.
- Record numbers of strains or isolates to be examined on the corresponding strips.
- Insert the rest of the unused strips into a storage bag with desiccant and store it in the refrigerator. Protect the plate against humidity. We recommend to use the plate within 4 weeks from the opening.

Note:

Any uneven distribution of substrate in the well does not affect the functionality of the test.

Inoculation:

- Homogenise the bacterial suspension in physiological solution well.
- Inoculate 0.1 ml of the suspension into all wells of a strip.
- Add two drops of paraffin oil to the wells H, G, F, E a D of the first row (tests URE, ARG, ORN, LYS and AAM).

Note:

The lid is labelled by test abbreviations and a graphic symbol ● for the tests with added paraffin oil to facilitate the inoculation and the reading. If you use the lid to cover the plate during the work, clean its inner side with ethanol before the use.

Incubation:

- Insert the frame with inoculated strips into a polyethylene bag.
- Fold the open end of the bag under the plate to prevent evaporation during the incubation.
- Insert the plate of NEFERMtest 24 into an incubator set up for 37 °C. Incubate it for 24 hours, extend the incubation for further 24 hours for slowly growing strains.

Note: In case of *Pseudomonas putida* identification, extend the incubation to 48 hours as it can reveal a metabolically less active *Pseudomonas aeruginosa*.

Evaluation:

Read all the tests after incubation and record the results into a record sheet.

It is the best to use the Colour scale for NEFERMtest 24 for assessment of colour reactions. You can use the table **Interpretation of reactions** (see in the Instructions for use below) or colour reactions of control strains to read the results.

Identification:

- Use the identification program ErbaExpert or the Code Book for NEFERMtest 24 for identification.
- Take into the consideration all the aspects of a culture: a source of the isolate, colony morphology, microscopy etc. Other characteristics are important in the identification of Gram-negative non-fermenting bacteria: pigment formation and fluorescence, growth at 42°C, hemolysis, growth on MacConkey agar, catalase, motility, gelatin test, results of additional tests OXItest and OFtest, etc.
- If the identification is not successful, repeat the NEFERMtest 24, or perform recommended additional tests.

Note:

Record sheets are designed to create so called profile. It is a numerical code that enables to find a result of the identification in a Code book easily. The procedure how to create a profile is described in the Code book.

Disposal of the used material:

- Insert the used plate into a jar for infectious material. Then autoclave it or incinerate it.
- Discard used paper packaging into recycling.

The most frequent causes of identification failure:

- Mixed or contaminated culture.
- Using inoculum of low density or low volume.
- Inoculum has contaminated adjacent strips used for another tested strain.
- The corresponding tests were not overlaid by paraffine oil.
- Failure to follow some of the steps in the recommended procedure.
- There may be an atypical strain or species or a representative of a species or a related genus present that is not included in the **Identification table**.

Performance of the kit:

The kit was tested on a set of 104 strains.

- The identification of 82,6% of the strains was correct.

Test quality control:

Quality of chemicals used for production of NEFERMtest 24 plates is checked by standard testing procedure. A functionality check for manufactured batches of plates is carried out on bacterial control strains.

If you want to perform quality control testing at your workplace we recommend to use the following control strains:

- *Pseudomonas aeruginosa* CCM 1960 (ATCC 10145)
- *Myroides odoratus* CCM 3296 (ATCC 4651)
- *Burkholderia cepacia* CCM 3919
- *Enterobacter (Cronobacter) sakazakii* CCM 3461

These strains are supplied by the CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz. The strains are supplied lyophilised or on gelatin discs.

Caution:

Only fresh isolates of CCM strains must be used to check the functionality of the kit. **These CCM control strains were selected to check the functionality of the kit only. They should not be used to check the accuracy or success of the identification.**

Health protection: Components of the kit are not classified as dangerous.

Interpretation of reactions:

Column	Test	Code	Reaction	
			positive	negative
Row 1 (1st, 4th, 7th, 10th row of the plate)				
H	Urease	URE	red, red to orange	yellow, light orange
G	Arginin	ARG	purple, blue	green
F	Ornithin	ORN	blue, blue to green	yellow to green, green
E	Lysin	LYS	blue, blue to green	yellow to green, green
D	Acetamid	AAM	blue to green, green	yellow, yellow to green
C	β - Glukosidase	bGL	yellow, yellowish	colourless
B	N-acetyl-β-D-glucosaminidase	NAG	yellow, yellowish	colourless
A	Simmons citrate	SCI	blue, blue to green	yellow to green, green
Row 2 (2nd, 5th, 8th, 11th row of the plate)				
H	Laktose	LAC	yellow, yellow to grey	violet, violet to grey
G	Mannitol	MAN	yellow, yellow to grey	violet, violet to grey
F	Trehalose	TRE	yellow, yellow to grey	violet, violet to grey
E	Xylosa	XYL	yellow, yellow to grey	violet, violet to grey
D	Arabinose	ARA	yellow, yellow to grey	violet, violet to grey
C	β – Galaktosidase	aGA	yellow, yellowish	colourless
B	β – Galaktosidase	bGA	yellow, yellowish	colourless
A	Malonate	MAL	blue, blue to green	yellow, green
Row 3 (3th, 6th, 9th, 12th row of the plate)				
H	Galaktose	GAL	yellow, yellow to grey	violet, violet to grey
G	Maltose	MLT	yellow, yellow to grey	violet, violet to grey
F	Cellobiose	CEL	yellow, yellow to grey	violet, violet to grey
E	Saccharose	SUC	yellow, yellow to grey	violet, violet to grey
D	Inositol	INO	yellow, yellow to grey	violet, violet to grey
C	γ -Glutamyltransferase	gGT	yellow, yellowish	colourless
B	Phosphatase	PHS	yellow, yellowish	colourless, cream-coloured
A	Esculine	ESL	black, dark brown	beige, light brown

Control strains:

Row	H	G	F	E	D	C	B	A
Pseudomonas aeruginosa CCM 1960								
1	URE	ARG	ORN	LYS	AAM	bGL	NAG	SCI
	+	+	-	-	+	-	-	+
2	LAC	MAN	TRE	XYL	ARA	aGA	bGA	MAL
	-	d	d	+	+	-	-	+
3	GAL	MLT	CEL	SUC	INO	gGT	PHS	ESL
	+	-	-	-	-	+	s	-
Myroides odoratus CCM 3296								
1	URE	ARG	ORN	LYS	AAM	bGL	NAG	SCI
	+	-	-	-	-	-	+	-
2	LAC	MAN	TRE	XYL	ARA	aGA	bGA	MAL
	-	-	-	-	-	-	-	-
3	GAL	MLT	CEL	SUC	INO	gGT	PHS	ESL
	-	-	-	-	-	-	+	-
Burkholderia cepacia CCM 3919								
1	URE	ARG	ORN	LYS	AAM	bGL	NAG	SCI
	-	-	+	+	d	-	+	s
2	LAC	MAN	TRE	XYL	ARA	aGA	bGA	MAL
	+	d	+	d	d	d	d	d
3	GAL	MLT	CEL	SUC	INO	gGT	PHS	ESL
	+	+	+	-	d	+	+	-
Enterobacter (Cronobacter) sakazakii CCM 3461								
1	URE	ARG	ORN	LYS	AAM	bGL	NAG	SCI
	-	+	d	-	-	+	+	+
2	LAC	MAN	TRE	XYL	ARA	aGA	bGA	MAL
	+	+	+	+	+	+	+	-
3	GAL	MLT	CEL	SUC	INO	gGT	PHS	ESL
	+	+	+	+	d	+	s	+

Explanations: + = positive reaction - = negative reaction s = weak reaction d = variable reaction

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