Kat. č.: MLT00011

**Pro mikrobiologii**

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Cat. No.: MLT00011

**For microbiology**

The NEISSERIAtest is a miniaturized version of conventional procedures for the identification of *Neisseria* species. It is a ready-to-use microwell plate system designed for performance of 7 biochemical tests (acid production from glucose, maltose, fructose and sucrose; detection of  $\gamma$ -glutamyl transferase, hydrolysis of tributyrin and synthesis of polysaccharide). Used together with some supplementary tests (e.g. cytochrome oxidase,  $\beta$ -galactosidase) enable identification of 10 species of *Neisseria* and *Moraxella (Branhamella) catarrhalis* accurately and easily.

- The kit NEISSERIAtest contains:**
- 3 microtitration plates (for the identification of 12 strains each) with desiccant
  - Instructions for use including the differentiation table,
  - Colour scale for NEISSERIAtest
  - 3 polyethylene bags for incubation
  - Storage bag (for storage of an open plate), 1 pc
  - 36 record sheets
  - Lid

**Storage, expiration:**

The NEISSERIAtest kit should be stored at (+2 to +8) °C. The expiration date is indicated on the package.

**Recommended procedure****Material required to perform a test (not included in the NEISSERIAtest kit):**

- Suspension medium for NEISSERIAtest (Cat. No. MLT00025 – 18 determinations)
- Petri dishes with cultivation medium suitable for neisseria
- Lugol iodine solution (Iodine powdered crystals 5 g; potassium iodine 10 g; distilled water 100 ml)
- Automatic micropipette (set to 0.1 ml) plus sterile tips
- Candle extinction jar or CO<sub>2</sub> incubator (4–10 % CO<sub>2</sub>) 37°C
- Inoculating loop, burner, marking pen
- Instrument DENSILAMETER II (Cat. No. INS00062)

**Material required to perform a supplementary test (not included in the NEISSERIAtest kit):**

- OXltest strips (Cat. No. MLT00039 – 50 determinations)
- Reagent for OXIDASE test (Cat. No. MLT00022 – 250 determinations)
- ONPtest strips (Cat. No. MLT00038 – 50 determinations)
- Sterile test tubes (100 x 15) mm with 1 ml of saline for ONPtest
- V+K DISK (Cat. No. MLT00087 - 100 determinations)

**For results evaluation (not included in the kit):**

- Code Book for NEISSERIAtest - located at [www.erbalachema.com](http://www.erbalachema.com)
- The ErbaExpert Identification Program

**Caution:**

For professional use only

**Follow the principles for working with infectious material!**

**Isolation of cultures:**

Perform the isolation of culture by usual techniques on recommended media (e.g. Mueller-Hinton chocolate agar, Thayer-Martin agar).

Before preparing a suspension of the culture in suspension medium the following should be checked:

- Cell morphology and Gram stain reaction
- Catalase and cytochrome oxidase tests

The culture should comprise Gram negative diplococci with positive catalase and cytochrome oxidase tests.

**Note:**

For the selective isolation of *Neisseria meningitidis* from the throat diagnostic disc V+K DISK (vancomycin + colistin), Cat. No. MLT00087 can be used. In the zone of antibiotic diffusion grow only *N. meningitidis*, *N. lactamica* and yeasts. Their differentiation:

- yeasts – by microscopy
- *N. lactamica* – by ONPtest, Cat. No. MLT00038 – up to 1 h ONPG gives positive reaction (yellow colour), ONPG negative strains should be identified by NEISSERIAtest
- N. meningitidis* – (ONPG negative) – to be identified by NEISSERIAtest with the result up to 4 h

**Preparation of the NEISSERIAtest plate:**

- Open the aluminium sachet close to the weld and take out the plate.
- Cut off required number of triple strips from the plate.
- Remove the adhesive tape from individual strips and insert them into prepared frame. In case you work with MIKROLATEST® kit for the first time and an empty frame is not available, use the frame of the first plate. The unused strips of the first plate put into the storage bag freely.
- Record number of the strains or isolates to be examined on the appropriate strips.
- Put the rest of the plate with desiccant to the bag for storage of open plate, enclosed with the kit, and store it in a refrigerator for further use; keep in mind to protect it from humidity. It is recommended to spend the rest of plate till 4 weeks after first use.

**Note:** Any uneven distribution of substrate in the well does not affect the functionality of the test.

**Preparation of inoculum:**

- Pre-warm the suspension medium and the NEISSERIAtest kit to (35–37) °C
- Use a pure, 24 hours old culture of Gram-negative, catalase and cytochrome oxidase positive diplococci.

The use of an older culture may lead to false reactions in some tests

- Open the ampoule with suspending medium by snapping off the top
- From the pure culture on Mueller-Hinton chocolate agar or Thayer-Martin agar make a suspension in the suspension medium having a turbidity equal to McFarland No. 3 turbidity scale
- Homogenize suspension thoroughly

Streak-out a sample from the inoculated suspension medium on a culture plate and incubate at (35–37) °C for 24 hours in a 4–10 % CO<sub>2</sub> atmosphere. This will provide a test to confirm purity of the suspension used to inoculate NEISSERIAtest and ONPtest.

**Inoculation:**

- Inoculate the strips as soon as possible; not later than 30 min after preparation of the suspension
- Homogenize the suspension thoroughly before use
- Inoculate 0.1 ml of the suspension into each well of the respective strip
- Do not dispense the inoculum too vigorously to prevent aerosol contamination of the wells of adjacent strip
- Cover the inoculated strips with lid
- Inoculation of the ONPtest: Insert ONPtest strip with the sterile forceps into tube with remaining inoculum or in inoculum made in saline

**Incubation:**

- Insert the frame with strips into a polythene bag
- Fold the open end of the bag under the plate to prevent desiccation during incubation
- Incubate NEISSERIAtest and ONPtest at (35–37) °C for 24 hours in aerobic conditions
- Preliminary reading of the results can be carried out already after 4 hours incubation
- The SPS test (well A) can be read only after 24 hours

**Reading:**

- Examine the purity of culture on control plate medium. Plates with no growth should be reincubated and examined after another 24 hours incubation.
- Record the results of all reactions (wells G–B) and the ONPtest (see note below).
- Compare the results of sugars (wells G–D) with negative control (well H). Any change of the colour reactions of these sugars in comparison to negative control in well H means positive reaction.
- Add one drop of the Lugol reagent into well A (test SPS) and record the result of the reaction.

**Note:**

- To evaluate the colour reactions of the tests follow the table "Interpretation of reactions" and/or the colour reaction of the control strains

**Identification:**

- For the identification refer to the "Differentiation table".
- If you have failed to identify the culture repeat the procedure as above.
- To complete the identification take into consideration all the results available, i. e. the source of culture, appearance and consistency of colonies, growth on selective media etc.

**The most frequent causes of identification failure:**

- Contaminated culture
- Using inoculum of low density or small volume of inoculum
- Inoculum has contaminated adjacent strip
- Failure to follow the recommended procedure
- There may be a species or strain whose data are not included in the "Differentiation table"

**Performance:**

- The kit was tested on a set of 112 strains.
- The identification of 98 % strains was correct.
- 2 % of the strains were not differentiated.

**Limitation:**

NEISSERIAtest should be used for the identification of *Neisseria* strains. Other bacteria do not necessarily give reliable results when using NEISSERIAtest.

**Disposal of used material:**

- After use, all ampoules, tips and strips must be autoclaved or incinerated
- Put paper packaging waste to recycling

**Quality control of NEISSERIAtest:** The quality control of the kits is performed systematically at various stages of their production. The batches are checked by tests on standard bacterial cultures. For those who wish to perform their own quality control tests, cultures mentioned in the table **Control strains** are recommended. **These strains are used to check the functionality of the kit, not to check the accuracy or success of the identification!**

*Neisseria subflava* CCM 3482;

*Neisseria lactamica* CCM 4392;

*Moraxella (Branhamella) catarrhalis* CCM 4391.

These strains are supplied in freeze-dried ampoules by the CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz)

## Control strains

| Strain                            | CCM  | Well. No. |          |          |          |          |          |          |          | Strip ONPtest |
|-----------------------------------|------|-----------|----------|----------|----------|----------|----------|----------|----------|---------------|
|                                   |      | H<br>NEC  | G<br>GLU | F<br>MLT | E<br>FRU | D<br>SUC | C<br>GGT | B<br>TRB | A<br>SPS |               |
| <i>Neisseria subflava</i>         | 3482 | –         | +        | +        | +        | +        | +        | –        | +        | –             |
| <i>Neisseria lactamica</i>        | 4392 | –         | +        | +        | –        | –        | –        | –        | –        | +             |
| <i>Moraxella (B.) catarrhalis</i> | 4391 | –         | –        | –        | –        | –        | –        | +        | –        | –             |

**Explanations:** +.....positive reactions;  
 –.....negative reactions;  
 NEC....negative control

**Caution:** It is necessary to use fresh isolates of the CCM strains each time when kit checked up for functionality. These strains are intended for check-up of the functionality of the kit, not for check-up of accuracy or effect of the identification.

**References:** Morse, S. A., and Knapp, J.S.: The genus *Neisseria*. In: The Prokaryotes. 2<sup>nd</sup>. ed. Edited by E. Ballows, H. G. Trüpper, M. Dworkin, W. Harder, and K. H. Schleifer. Vol. 3. p. 2495–2529. Springer Verlag, New York, Berlin, Heidelberg 1991.

**Health protection:** Components of the kit do not contain any dangerous substances.

### Interpretation of reaction

| Well No. | Test                     | Code | Reaction                       |                               |
|----------|--------------------------|------|--------------------------------|-------------------------------|
|          |                          |      | Positive                       | Negative                      |
| H        | Negative control         | NEC  | –                              | red to orange-red             |
| G        | Glucose                  | GLU  | yellow or pale orange          | red to orange-red             |
| F        | Maltose                  | MLT  | yellow or pale orange          | red to orange-red             |
| E        | Fructose                 | FRU  | yellow or pale orange          | red to orange-red             |
| D        | Sucrose                  | SUC  | yellow or pale orange          | red to orange-red             |
| C        | λ-glutamyl transferase   | GGT  | yellow or pale yellow          | colourless                    |
| B        | Tributyrin               | TRB  | yellow or yellow-orange        | red or orange-red             |
| A        | Polysaccharide synthesis | SPS  | black, dark blue or dark brown | yellow, brown or brown-yellow |
| Strip    |                          |      |                                |                               |
| OXItest  | cytochrome oxidase       | OXI  | dark blue or blue              | grey                          |
| ONPtest  | β-galactosidase          | ONP  | yellow or pale yellow          | colourless                    |

### Differentiation table

| NEISSERIAtest              |     |     |     |     |     |     |     |     |     | Supplementary tests |                                |   |
|----------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|---------------------|--------------------------------|---|
| Well No.:                  | H   | G   | F   | E   | D   | C   | B   | A   |     |                     | Growth on Thayer-Martin medium |   |
| Test:                      | NEC | GLU | MLT | FRU | SUC | GGT | TRB | SPS | ONP | DNA                 |                                |   |
| <i>N. gonorrhoeae</i>      | –   | +   | –   | –   | –   | –   | –   | –   | –   | –                   | –                              | + |
| <i>N. meningitidis</i>     | –   | +   | +   | –   | –   | +   | –   | –   | –   | –                   | –                              | + |
| <i>N. lactamica</i>        | –   | +   | +   | –   | –   | –   | –   | –   | +   | –                   | –                              | + |
| <i>N. polysaccharea</i>    | –   | +   | +   | –   | d   | –   | –   | +   | –   | –                   | –                              | + |
| <i>N. sicca</i> +          | –   | +   | +   | +   | +   | d   | –   | +   | –   | –                   | –                              | – |
| <i>N. mucosa</i> +         |     |     |     |     |     |     |     |     |     |                     |                                |   |
| <i>N. subflava</i>         | –   | +   | +   | d   | d   | d   | –   | d   | –   | –                   | –                              | – |
| <i>N. flavescens</i>       | –   | –   | –   | –   | –   | –   | –   | +   | –   | –                   | –                              | – |
| <i>N. cinerea</i> ++       | –   | –   | –   | –   | –   | –   | –   | –   | –   | –                   | –                              | – |
| <i>N. elongata</i> ++      |     |     |     |     |     |     |     |     |     |                     |                                |   |
| <i>M. (B.) catarrhalis</i> | –   | –   | –   | –   | –   | –   | +   | –   | –   | +                   | –                              | – |

**Explanations:** + = 75–100 % of positive reactions  
 – = 0–25 % of positive reactions  
 d = 26–74 % of positive reactions  
 NEC = negative control  
 DNA = deoxyribonuclease

**Notes:** + *N. sicca* and *N. mucosa* are not differentiated by NEISSERIAtest  
 ++ *N. cinerea* cells occur as plump cocci arranged in pairs or in scattered clusters  
*N. elongata* cells occur as short and slender rods often arranged as diplobacilli or in short chains

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