Kat. č.: MLT00012

Pro mikrobiologii

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Cat. No.: MLT00012

For microbiology

The STAPHYtest 16 kit is intended for the differentiation of important species of the genus *Staphylococcus* and other grampositive catalase positive coccus within 24 hours. It is a ready-to-use microwell plate system with 16 biochemical tests. Kit enables sixty examinations.

The identification can be supplemented by the paper strip tests for the Voges-Proskauer reaction, cytochrome oxidase and pyrrolidonylarylamidase in the form of VPtest, OXItest and PYRAtest. These tests are provided separately.

The STAPHYtest 16 kit contains:

- 10 microtitration plates (for identification of 6 strains each) with desiccant
- Instructions for use including the differentiation table
- Colour scale for STAPHYtest 16 kit
- 10 polyethylene bags for incubation
- Storage bag (for storage of open plate), 1 pc
- 60 record sheets
- Lid

Storage, expiration:

The STAPHYtest 16 kits should be stored in a refrigerator at (+2 to +8) °C. Expiration is indicated on the package.

Recommended working procedure for STAPHYtest 16

**Material required for work
(not included in the kit):**

- Reagent for PHOSPHATASE (Cat. No. MLT00018 – 250 determinations)
- Reagent for NITRATE (Cat. No. MLT00021 – 460 determinations)
- Paraffin oil, sterilized (Cat.No. MLT00042 – 250 determinations)
- VPtest (Cat. No. MLT00041 – 50 determinations)
- Reagent for ACETOINtest (Cat. No. MLT00016 – 90 determinations)
- Petri dishes with blood agar
- Sterile test tubes (100x 15) mm with 2.6 ml of sterile physiological saline
- Instrument DENSILAMETER II (Cat. No. INS00062)
- Vortex V1 (Cat. No. 50001715)
- Pipette 0.1 ml, sterile tips
- Thermostat 37 °C
- Usual microbiological laboratory equipment (loops, markers, burner)

**Material required for work with additional tests
(not included in the kit):**

- OXItest (Cat. No. MLT00039 – 50 determinations)
- Reagent for OXIDASEtest (Cat. No. MLT00022 – 250 determinations)
- PYRAtest (Cat. No. MLT00040 – 50 determinations)
- Reagent for PYRtest (Cat. No. MLT00023 – 800 determinations)

**Evaluation equipment
(not included in the kit):**

- Code Book for STAPHYtest 16 - located at www.erbalachema.com
The ErbaExpert Identification Program

Note:

If Code book is used for evaluation it is necessary to perform VPtest to create correct profile.

Caution:

- For professional use only

Follow principles for working with infectious material!

Isolation of cultures:

- Perform the isolation of culture by usual technique.
- From a pure culture perform Gram staining and catalase test.
- Use Gram-positive, catalase-positive cocci for identification by STAPHYtest 16.

Preparation of inoculum:

- From a pure culture on blood agar make a suspension in saline.
- The suspension must have a turbidity equal to McFarland No. 2 turbidity scale.

Culture purity control:

If required confirm purity of the suspension by streaming-out a sample from the inoculated suspension medium on cultivation medium. Incubate at 37 °C. Check after 24 hours, event. after 48 hours.

- Preparation of STAPHYtest 16 plate:**
- Open the aluminium sachet close to the weld and take out the plate.
 - Cut off required number of strips from plate.
 - Remove the adhesive tape from individual strips and insert them into prepared frame. In case, that you work with MIKROLATEST® kit for the first time and an empty frame is not available, use the frame of the first plate. The unused strips of the first plate put into the storage bag freely.
 - Record number of the strains or isolates to be examined on the appropriate strips.
 - Put the rest of the plate with desiccant to the storage bag, enclosed with the kit, and store it in a refrigerator for further use; keep in mind to protect it from humidity. It is recommended to spend the rest of plate till 4 weeks after first use.

Note:

Any uneven distribution of substrate in the well does not affect the functionality of the test.

Inoculation:

- Homogenize suspension thoroughly.
- Inoculate 0.1 ml of the suspension into wells.
- After inoculation, overlay the following tests with paraffin oil: microwells H, G and F (row 1, tests urease, arginine, ornithine).
- Inoculation of VPtest: Insert VPtest strip with the sterile forceps into tube with remaining inoculum (volume 1.0 ml).
- Furthermore is recommended to make disc test of sensitivity of bacitracin, novobiocin and furazolidon for differentiation of staphylococcus and micrococcus.

Note: a lid is printed with abbreviated names of tests and graphic symbols:

- (drop paraffine oil) and ◊ (drop reagent).
- Clean the inside of the lid by ethanol just before usage.

Note:

- For checking the correct course of tests and for comparing of colour reactions in the reading results the inoculation of control strains is recommended.
- For the performance of VPtest and OXItest follow the instructions in appropriate leaflets.

Incubation:

- Put the frame with inoculated strips into a polythene bag.
- Fold the open end of the bag under the plate to prevent dessication during incubation.
- Incubate the inoculated control plate and STAPHYtest 16 at 35–37°C for 24 hours; incubate VPtest strip at 35–37°C for 1.5 hour only.

Reading:

- After 1.5 hours of incubation read the VPtest result:
 - Add 3 drops of reagent VPT I and 3 drops of reagent VPT II.
 - Shake the tube thoroughly and incubate for another 30 min. at 35–37°C and record the results.
 - After 24 hours of incubation add one drop of reagents into the following wells:
 - well B (NITRATE) – reagent NIT
 - well A (PHOSPHATASE) – reagent PHS
- Record the results of all reactions. Read the reactions according to Colour Scale for STAPHYtest 16 or use table “Interpretation of reactions” in this instruction or compare with the colour reactions of the control strains.

Note:

- Into wells B with negative test for nitrate add carefully small amount of zinc powder (about 0.5 mg) to confirm negative reaction; in case of negative reaction red colour will appear within 10 min.
- Read the tests bGA and GLR (row 1; wells E and D) against white background.
- For confirmation and increase of identification it is recommended to make a test for detection of activity of pyrrolidonylarylamidase – PYRAtest: follow included working information.
- For the identification of veterinary strains (for group of *S. sciuri* / *S. lentus*, ev. *M. caseolyticus*), it is possible to use OXItest which is provided separately. Perform OXItest according to instruction for OXItest.

Identification:

- For the identification use the identification software ErbaExpert, Code book or “Identification table” in this instruction.
- To complete the identification take into consideration all the results including additional characteristics available, i. e. the source of culture, appearance and consistency of colonies etc.
- If you have failed to identify the culture repeat the procedure as above.

The most frequent causes of identification failure:

- Contaminated culture.
- Using inoculum of low density or small volume.
- Inoculum has contaminated adjacent strip.
- The corresponding tests were not overlaid by paraffin oil.
- The reagent was dropped into adjacent row.
- Failure to follow the recommended procedure.
- There may be a species or strains whose data are not included in the “Identification table” or Register.

Disposal of used material:

- After use, all ampoules, strips and tips must be autoclaved or incinerated.
- Put paper packaging waste to recycling.

Performance:

- The kit was tested on a set of 117 strains.
- The identification of 94.9% strains was correct.
- The identification of 3.4% strains was correct, but with low % i.d.
- 1.7% of the strains were not identified.

Quality control of STAPHYtest 16:

The kits are systematically quality controlled at various stages of their manufacture. The batches are checked by means of standard bacterial cultures. For those who wish to perform their own quality control tests, cultures mentioned in the table below are recommended. These strains are intended for check-up of the functionality of the kit, not for check-up of accuracy or effect of the identification!

Test results of control strains

CCM No.	Row 1								Row 2							
	H URE	G ARG	F ORN	E bGA	D GLR	C ESL	B NIT	A PHS	H GAL	G SUC	F TRE	E MAN	D XYL	C MLT	B MNS	A LAC
2737	-	+	-	-	+	-	+	-	s	+	+	-	-	+	-	-
4069	-	-	+	-	-	-	+	-	-	+	+	-	-	+	+	-
4296	+	-	-	s	+	-	-	-	-	-	+	+	-	+	+	+
3572	+	-	-	+	+	+	+	+	+	+	s	+	+	+	+	+

Explanations: + = positive reaction - = negative reaction s = weak reaction

- *Staphylococcus haemolyticus* CCM 2737 (ATCC 29970)
- *Staphylococcus lugdunensis* CCM 4069
- *Staphylococcus cohnii* subsp. *urealyticum* CCM 4296 (ATCC 49331)
- *Staphylococcus gallinarum* CCM 3572 (ATCC 35539)

These strains are supplied in freeze-dried ampoules by the CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz

It is necessary to use fresh isolates of the CCM strains each time when kit checked up for functionality.

Interpretation of reactions

Column	Test	Code	Reaction	
			positive	negative
Row 1				
H	Urease	URE	red-violet, orange-red	yellow, pale orange
G	Arginine	ARG	red-violet, red	yellow, pale orange
F	Ornithine	ORN	red-violet, red	yellow, pale orange
E	β-Galactosidase	bGA	yellow, pale yellow	colourless
D	β-Glucuronidase	GLR	yellow, pale yellow	colourless
C	Esculin	ESL	black, dark brown	colourless, pale brown
B	Nitrate	NIT	dark red, red	colourless, pale pink
A	Phosphatase	PHS	red-violet	colourless, pale pink
Row 2				
H	Galactose	GAL	yellow, yellow-brown	violet, grey-violet
G	Sacharose	SUC	yellow, yellow-brown	violet, grey-violet
F	Trehalose	TRE	yellow, yellow-brown	violet, grey-violet
E	Mannitol	MAN	yellow, yellow-brown	violet, grey-violet
D	Xylose	XYL	yellow, yellow-brown	violet, grey-violet
C	Maltose	MLT	yellow, yellow-brown	violet, grey-violet
B	Mannose	MNS	yellow, yellow-brown	violet, grey-violet
A	Lactose	LAC	yellow, yellow-brown	violet, grey-violet
VPtest	Acetoin	VPT	red, pink	colourless, pale pink
OXItest	Oxidase	OXI	blue	colourless
PYRAtest	Pyrrrolidonylarylamidase	PYR	red, orange	yellow

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Explanation: + = 90–100% of positive reactions
 (+) = 75–89% of positive reactions
 d = 26–74% of positive reactions
 (-) = 11–25% of positive reactions
 - = 0–10% of positive reactions

1) = The taxons can be found in human clinical material.
 2) = The taxons can be rarely found in human clinical material.

Health protection: Components of the kit do not contain dangerous substances.

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