STREPTOtest 16

Kat. č.: MLT00014

Pro mikrobiologii

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Cat. No.: MLT00014

For microbiology

The kit STREPTOtest 16 is intended for the routine identification of streptococci, especially of the species in the humane clinical material. It is a ready-to-use microwell plate system with 16 biochemical tests. Kit enables sixty examinations.

The identification is complemented by tests supplied in form of detection strips: PYRAtest for the detection of the pyrrolidonylarylamidase activity (PYRase, test PYR) and VPtest for the detection of the acetoin formation.

The kit STREPTOtest 16 contains:

- 10 microtitration plates (for identification of 6 strains each) with desiccant.
- Instructions for use including the differentiation table
- Colour scale for STREPTOtest 16
- 10 polyethylene bags for incubation
- Storage bag (for storage of open plate), 1 pc
- 60 record sheets
- Lid

Storage, expiration:

STREPTOtest 16 should be stored at (+2 to +8) °C. Expiration is indicated on each package.

Recommended working procedure for STREPTOtest 16

Material required for work (not included in the kit):

- Suspension medium for STREPTO test 16 (Cat. No. MLT00026 – 20 determinations)
- Paraffin oil sterilized (Cat. No. MLT00042 – 370 determinations)
- Reagent for HIPPURATE test (Cat. No. MLT00019 – 400 determinations)
- Reagent for PHOSPHATASE test (Cat. No. ML00018 – 250 determinations)
- Petri dishes with the blood agar
- Test tubes (100x15) mm with 3 ml of a sterile physiological solution
- Instrument DENSILAMETER II (Cat. No. INS00062)
- Test tubes rack
- Automatic micropipette 0.1 ml, sterile tips
- Thermostat (35–37) °C
- Equipment of a usual microbiological laboratory (loops, markers, burner)

Material required for work with additional tests

(not included in the kit):

VPtest (Cat. No. MLT00041 – 50 determinations)
 Reagent for ACETOIN test (Cat. No. MLT00016)
 PYRAtest (Cat. No. MLT00040 – 50 determinations)
 Reagent for PYR test (Cat. No. MLT00023)

Evaluation equipment (not included in the kit):

- Code Book for STREPTOtest 16 - located at www.erbalachema.com
- The Erba Expert Identification Program

Note:

If Code book is used for evaluation it is necessary to perform VPtest to create correct profile.

Caution:

- For professional use only

Follow principles for working with infectious material!

Isolation of cultures:

- Perform the isolation of individual cultures by using usual bacteriological technique on the blood agar (Columbia agar, blood agar no. 4).
- Evaluate the morphology of the pure culture, the haemolytical activity on the blood agar, perform the test for the detection of the pyrrolidonylarylamidase activity (test PYR, detection strip PYRAtest), resp. other tests confirming the presence of a group antigene (A, B, C, D, F, G) and record the results obtained in the record sheet.

Notes:

- Incubate the viridans streptococci and pneumococci in 5% CO₂ atmosphere.
- To identify enterococci, kit EN-COCCUStest is recommended.
- For β-haemolytic streptococci (streptococci with large colonies), the presence of a group antigen is the primary information whereas for pneumococci, it is the sensitivity to optochin and the solubility in bile.

Preparation of inoculum:

- Prepare a suspension of a pure 24 hours old culture in a physiological saline. Homogenize the suspension thoroughly.
- The suspension turbidity has to be equal to the McFarland's No. 3 turbidity scale.

Note:

- To obtain the necessary suspension density, it is recommended to use swabs for the suspension preparation.

Culture purity control:

If required confirm purity of the suspension by streaming-out sample from the inoculated suspension medium on cultivation medium. Incubate at 37 °C. Check after 24 hours

Preparation of the STREPTOtest 16 plate:

- Open the aluminium sachet close to the weld and take out the plate.
- Cut off required number of strips from plate.
- Remove the adhesive tape from individual strips and insert them into prepared frame. In case, that you work with MIKROLATEST® kit for the first time and an empty frame is not available, use the frame of the first plate. The unused strips of the first plate put into the storage bag freely.
- Record number of the strains or isolates to be examined on the appropriate strips
- Put the rest of the plate with desiccant to the bag for storage of open plate, enclosed with the kit, and store it in a refrigerator for further use; keep in mind to protect it from humidity. It is recommended to spend the rest of plate till 4 weeks after first use.

Note:

Any uneven distribution of substrate in the well does not affect the functionality of the test.

Inoculation:

- Homogenize the bacterial suspension thoroughly before use.
- Inoculate 0.1 ml of the suspension into the first 8 wells (the columns H-A of the first row, i.e. the tests HIP through URE).
- Pipette 1 ml of the suspension to the volume of one vial with the suspension medium for STREPTOtest 16 (1 ml).
- Homogenize the diluted suspension thoroughly.
- Inoculate 0.1 ml of the diluted suspension into the remaining eight wells (columns H-A of the second row, the tests MAN through RIB).
- Avoid the contamination of adjacent microwells during the inoculation process.
- After the inoculation, add 2 drops of the sterile paraffin oil into the wells of the column B and A in the first row (the tests ARG and URE).
- Insert the VPtest strip in the test tube containing the remaining suspension (1 ml of suspension), so that both strip zones were submerged in the suspension and provide the test tube with a stopper.

Note: a lid is printed with abbreviated names of tests and graphic symbols:

- (drop paraffine oil) and ◊ (drop reagent).

Clean the inside of the lid by ethanole just before usage.

Note:

- Inoculate control strains simultaneously with any series of unknown strains and when using a new batch of STREPTOtest 16 to control the colour expression of positive and negative reactions.

Incubation:

- Put the frame with inoculated strips into the polyethylene bag for incubation.
- Fold the open end of the bag under the plate to avoid desiccation.
- Incubate the VPtest for 2 hours and the STREPTOtest 16 plates for 24 hours at 35–37 °C.

Results, evaluation:

- Evaluation of VPtest: add always 3 drops of the reagent for VPT I and 3 drops for VPT II into the test tube, shake thoroughly and incubate for other 30 minutes at (35–37)°C. Then read the VP reaction and record the result.
- Evaluation of STREPTOtest 16:
 - drop the reagents into the following wells:
 - 1st row, well H (hippurate test) – 2 drops of the reagent for HIP
 - 1st row, well G (phosphatase test) – 1 drop of the reagent for PHS.
 - let the plate incubate for 5–10 minutes at the laboratory temperature for correct development of colour reaction.
 - read all tests and record the results into the record sheet.

Note:

- Use the Colour scale for STREPTOtest 16, the table "Interpretation of reactions" or proceed in accordance with colour reactions of the control strains to evaluate the colour reactions properly.
- It is recommended to read the tests LAP, GLR and aGA against white background.
- The test HIP should be interpreted within 10 minutes after the reagent addition, after this period, false positive reactions might appear.

Identification:

- Use the Code book or Identification software ErbaExpert.
- To complete the identification evaluate the culture as a complex, considering the serological group, the haemolysis, the kind of colonies, the microscopy, the origin of the isolate as well as other typical characters.
- In the case of an identification failure, repeat STREPTOtest 16 and complement the identification by other tests.

Disposal of the used material:

- Put the used plate into the vessel intended for the infectious material and autoclave or destroy it by incineration.
- Put paper packaging waste to recycling.

The most frequent causes of identification failure:

- A mixed or contaminated culture.
- Inoculum applied of a low density or a small volume.
- Inoculum was spread even over the adjacent row.
- The arginine and urease tests have not been overlaid by the paraffin oil.
- The reagent was dropped into the adjacent well.
- Failure to follow the recommended working procedure.
- It might be an atypical strain or representative of a species or an akin genus that has not been included in the Identification table.

Performance:

The kit was tested on a set of 60 clinically important strains. The identification of more than 90 % strains was correct.

Quality control of tests:

The kits are systematically quality controlled at various stages of their manufacture. The batches are checked by means of standard bacterial cultures. For those who wish to perform their own quality control tests, cultures mentioned in the table below are recommended. **These strains are intended for check-up of the functionality of the kit, not for check-up of accuracy or effect of the identification!**

CCM No.	Row 1								Row 2							
	H HIP	G PHS	F LAP	E GLR	D aGA	C ESL	B ARG	A URE	H MAN	G SOR	F TRE	E LAC	D RAF	C INU	B MLB	A RIB
4043	-	+	+	-	-	-	+	-	-	-	+	-	-	-	-	-
4617	+	-	+	+	-	+	s	-	+	+	+	+	-	s	-	s
1911	+	-	-	-	+	s	-	-	+	-	-	-	+	-	s	s

Legend:

+ = positive reaction - = negative reaction s = slightly positive reaction

- *Streptococcus constellatus subsp. constellatus* CCM 4043
- *Streptococcus uberis* CCM 4617
- *Aerococcus viridans* CCM 1911 (ATCC 29503)

These strains are supplied by the CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz
The strains are supplied in a lyophilized form.

Note:

For the purpose of checking the positive reaction of the URE test, it is possible to use the strain *Staphylococcus cohnii subsp. urealyticum* CCM 4296

Health protection:

Components of the kit are not classified as dangerous.

INTERPRETATION OF REACTIONS

Column	Test	Code	Reaction	
			positive	negative
Row 1				
H	Hippurate	HIP	blue	colourless, pale bluish
G	Phosphatase	PHS	red-to-violet, red	colourless, pale pinky
F	Leucine aminopeptidase	LAP	yellow	colourless
E	β-Glucuronidase	GLR	yellow	colourless
D	α-Galactosidase	aGA	yellow	colourless
C	Esculin	ESL	black, dark brown	colourless, pale brown
B	Arginine	ARG	red-to-violet, red	yellow, pale orange
A	Urease	URE	red, red-to-orange	yellow, pale orange
Row 2				
H	Mannitol	MAN	yellow, yellow-to-orange	red, orange-to-red
G	Sorbitol	SOR	yellow, yellow-to-orange	red, orange-to-red
F	Trehalose	TRE	yellow, yellow-to-orange	red, orange-to-red
E	Laktose	LAC	yellow, yellow-to-orange	red, orange-to-red
D	Raffinose	RAF	yellow, yellow-to-orange	red, orange-to-red
C	Inulin	INU	yellow, yellow-to-orange	red, orange-to-red
B	Melibiose	MLB	yellow, yellow-to-orange	red, orange-to-red
A	Ribose	RIB	yellow, yellow-to-orange	red, orange-to-red

STREPTOtest 16

IDENTIFICATION TABLE

Serol. group	STREPTOtest 16																		Identification
	Row 1									Row 2									
	H	G	F	E	D	C	B	A	A	H	G	F	E	D	C	B	A	A	
EM	YR	PT	IP	HS	PL	GA	ER	AL	MS	OR	TE	LA	RA	IN	ML	RI	MB		
A	+	+	-	-	+	+	-	-	+	-	-	+	+	-	-	-	-	S. pyogenes	
B	(+)	-	+	+	+	d	-	-	+	-	-	(+)	d	-	-	-	(+)	S. agalactiae	
C	-	-	-	-	+	+	+	-	-	+	-	-	d	+	+	-	-	S. dysgalactiae	
C	+	-	-	-	+	+	+	-	-	+	-	-	-	-	-	-	-	S. equi subsp. equi	
C	+	-	-	-	+	+	+	-	(-)	+	-	-	+	+	-	-	(+)	S. equi subsp. zooepidemicus	
C	+	-	-	-	+	+	(+)	-	(-)	+	-	-	+	d	-	-	+	S. equisimilis	
G	+	-	-	-	+	+	(-)	d	d	+	-	-	(-)	(+)	-	-	+	S. group G	
L	+	-	-	d	+	+	+	-	-	+	-	-	+	(+)	-	-	+	S. group L	
E,P,U,V	+	-	+	+	+	+	d	+	+	-	+	(+)	+	d	-	-	(+)	S. porcinus	
R,S,RS,T	-	d	-	-	+	+	(+)	(+)	+	-	-	-	+	+	d	d	-	S. suis	
-,F,C,(G,A)	d	-	+	-	+	+	-	d	+	+	-	(-)	-	+	(+)	d	(-)	S. anginosus	
-,F(C,G)	d	-	+	-	+	+	-	-	+	+	-	-	(+)	(-)	-	-	d	S. constellatus	
-,F,C(G)	d	-	+	-	+	+	-	-	+	+	-	-	+	+	-	-	-	S. intermedius	
	-	-	+	-	-	+	-	(+)	+	-	-	+	+	+	+	(+)	(+)	S. mutans "group"	
	-	-	+	-	d	+	-	-	+	-	d	-	-	d	d	(+)	d	S. salivarius	
	-	-	+	-	-	+	-	-	d	-	+	-	-	(+)	-	-	-	S. vestibularis	
D	-	-	+	-	-	+	-	-	d	+	-	-	+	+	+	(+)	d	S. bovis I	
D	-	-	+	-	-	+	-	-	+	+	-	-	-	(-)	+	+	d	S. bovis II/1	
D	-	-	+	-	-	+	d	(+)	+	-	-	-	-	+	+	d	(-)	S. bovis II/2	
D	-	-	+	-	-	+	-	-	+	-	-	(-)	-	(-)	d	(-)	(-)	S. equinus	
-, E	-	(-)	+	+	d	+	+	(-)	+	+	-	+	+	+	+	(-)	+	S. uberis	
-, E	-	(-)	+	d	d	+	(-)	(-)	+	+	-	+	+	+	(-)	(-)	-	S. parauberis	
	-	(-)	-	+	d	+	d	-	(-)	(-)	-	d	-	(+)	(+)	-	-	S. acidominimus	
	-	(-)	-	-	+	-	+	(+)	(-)	d	-	-	-	+	+	(+)	d	S. pneumoniae	
	-	-	-	-	d	+	-	d	-	+	-	-	(-)	-	+	d	-	S. mitis 2	
	-	-	-	-	d	+	-	d	-	-	-	-	-	(-)	+	d	-	S. oralis / mitis 1	
	-	-	-	-	+	-	d	+	+	-	-	(-)	(+)	+	d	(+)	d	S. sanguis esl+	
	-	-	-	-	+	-	+	+	+	-	-	-	(+)	+	-	(+)	-	S. sanguis esl-	
	-	-	-	-	+	+	-	d	+	+	-	-	+	(+)	d	(+)	d	S. gordonii	
	-	d	(-)	(+)	-	-	(-)	d	(+)	-	-	d	(-)	+	(+)	d	(-)	Aerococcus viridans	

Legend:

- + = positive reaction
- (+) = mostly positive reaction
- HEM = β haemolysis
- d = variable reaction
- (-) = mostly negative reaction
- PYR = pyrrolidonylamidase (PYRAtest)
- VPT = acetoin (VPtest)
- = negative reaction

Note: In the case of a non-explicit differentiation of species, it is possible to classify the viridans streptococci into one of the following groups:

- S. salivarius group = S. salivarius, S. vestibularis
- „S. milleri“ group = S. anginosus, S. constellatus, S. intermedius
- S. sanguis group = S. sanguis, S. gordonii
- S. mitis group = S. mitis, S. oralis

(RUOFF, K. L. 1995. Streptococcus, p. 299–307. In: P.R. Murray et al. (ed.), Manual of Clinical Microbiology, 6th ed. ASM, Washington, D.C.)

The facultatively anaerobic, catalase negative, gram-positive cocci:

Genus	PYR	LAP	Vancomycine	Bile-esculin	NaCl 6.5%	Satellitism
Streptococcus	-	+	C	-	-	-
Enterococcus	+	+	C ³⁾	+	+	-
Lactococcus	+	+	C	+	d	-
Vagococcus	+	+	C	+	+	-
Globicatella	+	-	C	-	+	-
Leuconostoc	-	-	R	d	d	-
Pediococcus	-	+	R	+	d	-
Aerococcus ¹⁾	+	-	C	d	+	-
Gemella	+	d	C	-	-	-
Lactobacillus	-	-	C ³⁾	d	d	-
Helcococcus	+	-	C	+	+	-
Abiotrophia ²⁾	+	+	C	-	-	+

Note:

a simplified scheme for differentiation of akin genera

1) A. urinae = PYR-, LAP+

2) Formerly indicated as satellite, nutritionally deficient, thiol dependent streptococci

3) Representatives of the Enterococcus and Lactobacillus genus exhibit the vancomycine resistance in 1 % of cases

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