## INDOLtest

Kat. č.: MLT00028 **Pro mikrobiologii**

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Cat. No.: MLT00028

**For bacteriology**

INDOLtest is designed for rapid determination of indole, degradation product of bacterial metabolism of tryptophane. It is designed for a quick presumptive determination of *E. coli*, screening differentiation of indole positive and indole negative bacterial genera and as a necessary additional test for some identification kits.

**Principle:** Principle of the test is hydrolysis of tryptophan to indole, pyruvate and ammonia. Indole reacts with p - dimethylaminocinnamaldehyde (DMACA), the substrate in the solution of the INDOL-test, and produces a blue-green compound. INDOLtest is the most sensitive routine method for detection of bacterial tryptophanase. The sensitivity is 3 µg indol/ml.

**A package of INDOLtest contains:**

- 10 ml of solution for 140 determinations
- Instructions for use

**Storage, expiration:**

INDOLtest is sensitive to light and should be stored at (+2 to +8) °C in the dark. Pay attention to closing the bottle immediately after the use.

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For professional use only. The solution contains 15.6% hydrochloric acid 35-38%.

Warning

**Hazard statement:**

- H315 Causes skin irritation.
- H319 Causes serious eye irritation.
- H335 May cause respiratory irritation.

**Precautionary statement:**

- P280 Wear protective gloves/protective clothing/eye protection.
- P302+P352 IF ON SKIN: Wash with plenty of water.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

The packaging should be used only by trained individuals skilled in work with infectious material and in safety disposal according to local laboratory directives.

**Material required:**

- Pipette with sterile tips, filter paper (e.g. Whatman N. 1)
- Inoculation loop, burner, marking pen etc.

**Recommended procedure:**

- Use a 18-24-hour culture from a culture medium containing tryptone (tryptophane), preferably blood agar. We recommend a minimum tryptone content of 7 g/l in the medium for correct implementation of the test. If there is any doubt about the suitability of the medium, use a verified indole-positive strain to check the medium.
- Add a drop of the solution INDOLtest on a piece of filter paper so that the paper is sufficiently wet.
- Using an inoculation loop smear few colonies of the tested culture onto the drop zone on the filter paper. The paper must not dry out before the bacterial colonies are applied.
- Incubate at ambient temperature for 2 to 5 minutes.
- Evaluate the colour reaction directly on the drop zone according to the table **Interpretation of the reaction.**

**Interpretation of the reaction:**

Reaction	Colour expression
positive	blue-green
negative	pink

**Waste disposal:**

Used filter paper should be treated as potentially infectious material and should be disposed according to local and national regulations for the safe handling of such materials.

**The most frequent causes of identification failure:**

- Mixed or contaminated culture. If more than one strain is inoculated in the same culture plate a situation might occur that an indole positive strain can influence an originally indole negative strain into a slightly positive reaction. It is necessary to culture the strains on separate plates and repeat the INDOLtest.
- There is not enough tryptophane in the culture medium. Verify positive and negative reaction with recommended control strains.
- The culture media contains glucose that inhibits production of indole.
- The filter paper is not sufficiently moistened with INDOLtest solution.
- Violation of some steps of recommended procedure.
- They may be a species or strains for which INDOLtest is not recommended.

**Caution: Pay attention that the INDOLtest is not interchanged with Reagent for INDOL test (Cat. No. 10003372) that is an auxiliary colour forming reagent in some kits of the MIKRO-LA-TEST product range.**

**Quality control of INDOLtest:**

It is recommended to use the following control strains for the internal quality control:

- CCM 4225 (ATCC 35218) *E. coli* (positive control)
- CCM 5852 (ATCC 13882) *K. pneumoniae* subsp. *pneumoniae* (negative control)

These strains are distributed by Czech Collection of Microorganisms,

Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ

Tel: +420 549 491 430, Fax: +420 549 498 289

<http://www.sci.muni.cz/ccm>, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz)

The strains are delivered in freeze-dried form or in gelatine discs.

**Caution:** It is necessary to use fresh isolates of CCM strains to check the function of INDOLtest.

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