## β-LACTAMtest

IVD

CZ

Kat. č.: MLT00034

Pro mikrobiologii

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Cat. No.: MLT00034

For microbiology

β-LACTAM test is intended for the rapid acidimetric detection of bacterial β-lactamase activity. This acidimetric method is recommended for use in testing β-lactamase production by Neisseria species, Haemophilus species and Staphylococcus species.

### Principle:

β-lactamase enzyme hydrolyses benzylpenicillin contained in the zone of detection strip, which is detected by colour change of pH indicator in the detection zone due to pH shift.

### A package of β-LACTAMtest contains:

- 50 test strips for 50 determinations
- Working leaflet

### Storage:

β-LACTAMtest should be stored in refrigerator at (+2 to + 8)°C. The expiration date is indicated on each package.

### Caution:

For professional use only. The packing should be used only by properly trained individuals being skilled in principles of work with infectious material and its safety disposals according to local obligatory directions of laboratory.

### Material required:

- Pipette, sterile tips
- Inoculation loop, burner, marking pen

### Recommended procedure:

- The isolation of cultures should be carried out by convection technique on Columbia blood agar (24 hour old culture)
- After taking the tube with strips out of the refrigerator let the tube before opening to equilibrate to ambient temperature to avoid forming condensate on the surface of tube and subsequently to protect the product from adverse effects of air humidity.
- Take out required amount of detection strips and immediately close the tube.
- It is recommended to spend the content of tube till 2 months after first use.
- Add condensation water from the lid of Petri dish on the reaction zone of the strip before smearing the bacterial mass on the strip. If no condensate water is present, moisten the reaction zone by 10 μl of water. Depending on the water pH value, it can happen, that the strip zone becomes red. This change has no influence on the identification.
- Using an inoculation loop smear several colonies of tested culture on the reaction zone of the strip. Also alternative procedure can be used: moisten the zone of strip and wipe off the colony directly from culture medium. Bacterial mass has to be spread but mustn't cover all coloration of the strip.
- Place the strip on suitable pad and incubate at ambient temperature for 10 minutes.
- Evaluate the colour reaction directly on the zone of the strip according to table **Interpretation of the reaction**.

**Note:** It is possible, that some strains of Staphylococcus aureus produce a small amounts of the β-lactamase enzyme unless enzyme production is induced by growth in a medium containing sub-inhibitory concentrations of β-lactam antibiotic. In this case the course of reaction can be significantly prolonged (till 60 min.).

### Interpretation of the reaction:

Reaction	Colour expression
positive	red
negative	yellow

### Waste disposal:

Strips should be treated as potentially infectious and should be disposed according to local and national regulations relating to the safe handling of such materials. After use strips must be autoclaved or incinerated. Packaging waste can be put to recycling or in municipal waste.

### The most frequent causes of identification failure:

- Contaminated culture.
- Violation of some steps of recommended procedure.
- There may be a species or strains for which β-LACTAMtest is not recommended.

### Standard control strains:

It is recommended to use standard control strains for checking of the reaction and interpretation of the colours of the reaction.

CCM 7221 <i>S. epidermidis</i>	Positive reaction
CCM 3953 <i>S. aureus</i>	Negative reaction

These strains are distributed by Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ Tel: +420 549 491 430, Fax: +420 549 498 289 <http://www.sci.muni.cz/ccm>, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz) The strains are delivered in a lyophilized form or on gelatinized disks.

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