# INDOXYLtest

(CZ)

Kat. č.: MLT00037

Pre mikrobiologii

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Cat. No.: MLT00037

For microbiology

INDOXYLtest is intended for the rapid detection of acetateesterase activity of *Moraxella catarrhalis* and other microorganisms.

**Principle:**

3 – indoxylacetate is hydrolyzed to leuko form of indigo. Free indoxyl groups reacts in presence of oxygen to chromogen form of indigo. Blue coloration of diagnostic zone represents positive reaction and demonstrates acetateesterase activity of examined microorganism.

**A package of INDOXYLtest contains:**

- 50 test strips for 50 determinations
- Working leaflet

**Storage, expiration:**

INDOXYLtest should be stored in refrigerator at (+2 to + 8)°C. The time of expiration is indicated on each package.

**Caution:**

For professional use only. The packaging should be used only by properly trained individuals being skilled in principles of work with infectious material and its safety disposals according to local obligatory directions of laboratory.

**Material required:**

- Pipette with sterile tips
- Inoculation loop, burner, marking pen etc.

**Recommended procedure:**

- The isolation of cultures should be carried out by convection technique on Columbia blood agar (24 hours old culture).
- After taking the tube with strips out of the refrigerator let the tube before opening to equilibrate to ambient temperature to avoid forming condensate on the surface of tube and subsequently to protect the product from adverse effects of air humidity.
- Take out required amount of detection strips and immediately close the tube.
- It is recommended to spend the content of tube till two months after first use.
- Add condensation water from the lid of Petri dish on the reaction zone of the strip before smearing the bacterial mass on the strip. If no condensate water is present, moisten the reaction zone by 10 µl of water.
- Using an inoculation loop smear several colonies of tested culture on the reaction zone of the strip.
- Place the strip on suitable pad and incubate at ambient temperature for 2 to 5 minutes.
- Evaluate the colour reaction directly on the zone of the strip according to table **Interpretation of the reaction**.

**Note:**

INDOXYLtest can be used also for differentiation of campylobacter species. In this case please incubate strip up to 30 min.

**Interpretation of the reaction:**

Reaction	Colour expression
Positive	Blue, blue - green
Negative	No colour change

**Waste disposal:**

Strips should be treated as potentially infectious and should be disposed according to local and national regulations to the safe handling of such materials. After use - strips must be autoclaved or incinerated. Packaging waste can be put to recycling or in municipal waste.

**The most frequent causes of identification failure:**

- Contaminated culture
- Violation of some steps of recommended procedure
- They may be a species or strains for which INDOXYLtest is not recommended

**Quality control of INDOXYLtest:**

It is recommended following control strain for internal quality control: CCM 4391 *Moraxella* ( B ) *catarrhalis* (positive reaction). Distilled water serve as negative control (instead of bacterial culture).

This strain is supplied by Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ Tel: +420 549 491 430, Fax: +420 549 498 289 http://www.sci.muni.cz/ccm, e-mail: ccm@sci.muni.cz  
The strains are delivered in freeze-dried form or on gelatine discs.

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