# OXItest

Kat. č.: MLT00039

## Pro mikrobiologii

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Cat. No.: MLT00039

## For microbiology

The OXItest is intended for the detection of bacterial cytochrome oxidase. A package of the OXItest contains 50 test strips enabling 50 determinations.

**Principle:**  
The activity of the cytochrome oxidase is detected by colour reaction of N, N-dimethyl-1,4-phenylenediamine with α-naphtol which form indophenol blue. The sensitivity of the reaction can be increased by Reagent for OXIDASE test.

**A package of the OXItest contains:**  
• 50 test strips for 50 determinations  
• working leaflet

**Storage:**  
In refrigerator at (+2 to +8) °C. The time of expiration is indicated on each package.

**Caution:**  
• For professional use only

**Recommended procedure:**  
Pick a single, well-isolated colony from the plating medium by an inoculation platinum or plastic loop (or full loop in case of pure culture) and spread the inoculum on the reaction zone of the strip. It is possible to print the reaction zone of the strip directly on the pure culture on the plating medium. Read the colour reaction immediately (within 1 minute); later false positive reactions are possible. Reagent for OXIDASE test (Cat. No. MLT00022) can be used for increasing of the sensitivity of the reaction: Add cca 10 µl of the reagent on the reaction zone of the strip before spreading the bacterial mass on the strip.

**Waste disposal:**  
Used strip should be treated as potentially infectious and should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging waste to recycling or in municipal waste.

### Interpretation of the reaction

Test	Code	Reaction	
β-galactosidase (ONPG)	OXI	positive	negative
		blue, pale blue	grey

It is recommended to use standard control strains for checking of the reaction and interpretation of the colours of the reaction.

### Standard control strains

Strain	CCM	ATCC	OXI
Enterobacter sakazaki	3461		-
Pseudomonas aeruginosa	1960	10145	+

These strains are distributed by  
Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ  
Tel: +420 549 491 430, Fax: +420 549 498 289 http://www.sci.muni.cz/ccm, e-mail: ccm@sci.muni.cz

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