Suspenní médium MIC G+

Kat. č.: MLT00071

Pro mikrobiologii

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Cat. No.: MLT00071

For microbiology

Suspension Medium MIC G+ is supplementary preparation for MIC G+ kit which is designed for antibiotic susceptibility testing. The Suspension Medium MIC G+ is supplied in glass tubes in the sterile conditions. The content of one tube (13 ml) is sufficient for preparation of the suspension of one strain (inoculation of one MIC G+ plate). The preparation of the suspension should be carried out according to the procedure recommended for MIC G+ test.

Content:

30 tubes with suspension medium MIC G+
Instructions for use

Storage, expiration:

The Suspension medium for MIC C+ should be stored at the temperature of (2 – 25)°C until the expiry date indicated on the package.

Quality control:

Prepare suspension of control strain with density of McF 0.5. Transfer 60 µl of suspension into 13 ml of suspension medium MIC G+ in tube. Incubate at 35±2°C for 16-20 hours.
Enterococcus faecalis CCM 4224 (ATCC 29212) – turbid growth

Warning:

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. No dangerous reagents are included in the kit.

Waste disposal:

The residue of the suspension medium should be disposed in accordance with local and national regulations related to the safe handling of such materials. Recycle packaging paper waste.

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