Kat. č.: MLT00089

Kat. č.: MLT00088

Kat. č.: MLT00090

X - FAKTOR
V - FAKTOR
X+V - FAKTOR

Pro mikrobiologii

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For microbiology

Diagnostic discs X - ФАКТОР, V - ФАКТОР и X+V - ФАКТОР are intended for simple routine differentiation of *Haemophilus* spp. from clinical material.

Principle:

Members of the genus *Haemophilus* are dependent on exogenous growth factors X (hemin) and V (NAD) in medium free of these factors. Strain to be tested after incubation shows growth only around the disc containing required essential factor.

A package contains:

- 100 diagnostic discs for 100 determinations
- working leaflet

Storage, expiration:

Discs should be stored in dark in a refrigerator at -18°C, until the expire date indicated on the package.

Required material:

- Agar culture media plates without blood (Mueller Hinton agar, Columbia Blood Agar Base)
- common laboratory equipment (loops, sterile forceps, burner)

Recommended procedure:

- Streak regularly fresh 24 h culture of strain to be tested on agar plate free of the growth factors
- Aseptically place disc with growth factors at a distance at least 2 cm of each other
- Incubate in an atmosphere containing (5-10)% CO₂ at 37°C for (18-24) h or 48 h if necessary.

Evaluation:

After overnight incubation are examined colonies around the discs:

- bacteria requiring growth factor X grow only around the discs X and X+V-ФАКТОР
- bacteria requiring growth factor V grow only around the discs V and X+V-ФАКТОР
- bacteria requiring both growth factors grow only around the discs X+V - ФАКТОР

Caution:

- For professional use only.
- After use all discs must be autoclaved or incinerated.

Differentiation table:

Strain	X-ФАКТОР	V-ФАКТОР	X+V-ФАКТОР	Haemolysis (blood agar)
<i>H. influenzae</i>	-	-	+	-
<i>H. parainfluenzae</i>	-	+	+	-
<i>H. haemolyticus</i>	-	-	+	+
<i>H. paraaerolyticus</i>	-	+	+	+
<i>H. aegyptius</i>	-	-	+	-
<i>H. ducreyi</i>	+	-	+	-
<i>H. haemoglobinophilus</i>	+	-	+	-
<i>H. segnis</i>	-	+	+	-
<i>H. paraphrophilus</i>	-	+	+	-
<i>H. aphrophilus</i>	-	-	-	-

Standard control strains:

Quality control of diagnostic discs X, V and X+V ФАКТОР may be provided with standard control strains:

Strain	Required factor
<i>Haemophilus influenzae</i> CCM 4457	X+V
<i>Haemophilus haemoglobinophilus</i> CCM 4648	X
<i>Haemophilus segnis</i> CCM 6052	V
<i>Haemophilus aphrophilus</i> CCM 3582	none

These strains are distributed by:

CCM - Czech Collection of Microorganisms, Masaryk University, Faculty of Science,

Kamenice 5, building A25, 625 00 Brno, CZ

Tel. +420 549 491 430, Fax +420 549 498 289, http://www.sci.muni.cz/ccm, e-mail: ccm@sci.muni.cz

Strains are supplied in freeze-dried form or on gelatin discs.

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